



Product description

Order number	73.666
Product description	Sample tube, suitable for Hitachi analyser, material: PS, 3 ml, (LxØ): 38 x 17 mm, transparent, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Hitachi analyser
Base shape	conical base

Size

Max. volume	3 ml
Diameter	17 mm
Diameter of opening	10 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Sample tube, suitable for Hitachi analyser,
transparent**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest
subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 5000
Piece(s) / pallet 150000
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.0017 kg
Weight of case 9.3 kg
EAN of inner box 4038917083225
EAN case 4038917002172

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT