



Product description

Order number	73.667
Product description	Sample tube, suitable for Abbott analyser AXSYM, material: PS, 2.5 ml, (LxØ): 53 x 12 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	Abbott analyser AXSYM
Base shape	conical base

Size

Max. volume	2.5 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	53 mm
Length of product	53 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 2500
Type of smallest bag
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 2500
Piece(s) / pallet 90000
Depth of case 398 mm
Width of case 398 mm
Height of case 265 mm
Case volume 0.042 cbm
Weight of product 0.0025 kg
Weight of case 6.78 kg
EAN of inner box 4038917083232
EAN case 4038917002189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

