



## Product description

Order number	10.252.027
Product description	Liquor/CSF-Monovette® , nominal volume: 10 ml, without preparation, collection and storage of diagnostic specimens, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap blue, with paper label, "Liquor/Cerebrospinal Fluid", sterile, 1 piece(s)/blister

## Product characteristics

Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	blue
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	collection and storage of diagnostic specimens
Example of use	CSF, collection, transport and storage of CSF/liquor
Base shape	round base
Volume of work	10 ml
Connection	Luer (m)

## Size

Sample volume	10 ml
Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Liquor/CSF-Monovette® , 10 ml, cap blue, (LxØ):**  
**102 x 15 mm, 1 piece(s)/blister**

**Page 2**

Piston plunger material                      Polystyrene (PS)

**Purity & certification**

Satisfies the requirement                      IATA, ADR  
Product category                                      In vitrodiagnostic device, CE  
CE certified    CE  
Purity standard                                        sterile  
Sterilisation    Electron irradiation  
Batched    yes

**Packaging**

Minimum order qty.                                500  
Type of smallest                                        blister  
subpackaging  
Piece(s) / inner box                                100  
Piece (s) / outer case                              500  
Piece(s) / pallet                                      13500  
Depth of box    351 mm  
Width of box    151 mm  
Height of box    171 mm  
Depth of case                                        778 mm  
Width of case                                        368 mm  
Height of case                                        195 mm  
Case volume    0.0558 cbm  
Weight of product                                  0.0077 kg  
Weight of case                                        5.46 kg  
EAN of inner box                                    4038917104487  
EAN case    4038917040334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**