



Product description

| | |
|---------------------|--|
| Order number | 63.504.027 |
| Product description | Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, PP, transparent, cap assembled, with plastic label, label/print: transparent/green, "CSF/Liquor", with graduations, cap assembled, natural, sterile, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------------|--|
| Label/ Print | with plastic label |
| Colour of print/label | transparent/green |
| Graduation | yes |
| Closure type | screw cap |
| Example of use | CSF, collection, transport and storage of CSF/liquor |
| Base shape | round base |
| Volume of work | 5 ml |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|----------------------|-------|
| Diameter | 13 mm |
| Length of product | 75 mm |
| Length including cap | 83 mm |
| Length excluding cap | 75 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|-------------------------------|------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 13500 |
| Depth of box | 351 mm |
| Width of box | 151 mm |
| Height of box | 171 mm |
| Depth of case | 778 mm |
| Width of case | 368 mm |
| Height of case | 195 mm |
| Case volume | 0.0558 cbm |
| Weight of product | 0.0039 kg |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 5 ml, round base, transparent

Page 3

| | |
|------------------|---------------|
| Weight of case | 3.56 kg |
| EAN of inner box | 4038917408561 |
| EAN case | 4038917408554 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT