



### **Product description**

Order number	95.813.001
Product description	Dispenser box, material: PS, lid assembled

### **Product characteristics**

Base shape	flat base
------------	-----------

### **Size**

Diameter of base	81 mm
Diameter of lid	83 mm
Diameter of opening	81 mm
Length of product	110 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white
Closure material	Polystyrene (PS)
Colour of cap	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      35  
Piece(s) / inner box                      35  
Piece (s) / outer case                      35  
Piece(s) / pallet                              1260  
Depth of case                                398 mm  
Width of case                                 398 mm  
Height of case                                300 mm  
Case volume                                 0.0475 cbm  
Weight of product                         0.0623 kg  
Weight of case                                2.77 kg  
EAN of inner box                            4038917121644  
EAN case                                        4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

