



### **Product description**

Order number	95.963.001
Product description	Wall holder, for Multi-Safe euroMatic® , material: stainless steel

### **Product characteristics**

Suitable for	Multi-Safe euroMatic®
--------------	-----------------------

### **Size**

Width of product	161 mm
Length of product	222 mm

### **Material & colours**

Product material	stainless steel
Colour of product	silver/white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      20  
Piece(s) / pallet                              2880  
Depth of case                              323 mm  
Width of case                                213 mm  
Height of case                                145 mm  
Case volume                                0.01 cbm  
Weight of product                        0.246 kg  
Weight of case                                5.24 kg  
EAN of inner box                        4038917233620  
EAN case                                      4038917233613

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

