




## Product Specification

Round tube holder for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)

Page 1



### Product description

Order number	95.963.006
Product description	Round tube holder, for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020), material: POM, (LxWxH): 130 x 40 x 88 mm

### Product characteristics

Suitable for	round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)
--------------	---

### Size

Width of product	40 mm
Height of product	88 mm
Length of product	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Round tube holder for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)

Page 2

### Material & colours

Product material Polyoxymethylene (POM)  
Colour of product silver/white

### Purity & certification

Product category no medical device | IVD  
Batched no

### Packaging

Minimum order qty. 1  
Type of smallest subpackaging case  
Piece(s) / inner box 1  
Piece (s) / outer case 10  
Piece(s) / pallet 1680  
Depth of case 288 mm  
Width of case 218 mm  
Height of case 150 mm  
Case volume 0.0094 cbm  
Weight of product 0.4003 kg  
Weight of case 4.23 kg  
EAN of inner box 4038917364300  
EAN case 4038917364294

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**