



Product Specification
Disposal container, Multi-Safe twin plus REC,
2,000 ml

Page 1

Product description

Order number 77.3884.220
Product description Disposal container, Multi-Safe twin plus REC, 2,000 ml, grey, with plastic label, label/print: yellow/red, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 161 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 80 piece(s)/carton

Product characteristics

Label/ Print with plastic label
Colour of print/label yellow/red
Example of use disposal of pointed and sharp objects
Absorbent liner yes
Shape of opening oval
Removal aid - turning aid Removal and rotation aid

Size

Max. volume 2000 ml
Diameter of base 112 mm
Height of product 161 mm

Material & colours

Product material recycled Polypropylene (recycled PP)
Colour of product grey
Closure material Polypropylene (PP)
Colour of cap red
Colour of lid red

Purity & certification

Satisfies the requirement TRBA 250, UN 3291, DIN EN ISO 23907
Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	80
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	80
Piece(s) / pallet	960
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.194 kg
Weight of case	17.17 kg
EAN of inner box	4066916004088
EAN case	4066916004095

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

