



Product Specification
Disposal container, Multi-Safe twin plus REC,
3,000 ml

Page 1

Product description

Order number	77.3884.230
Product description	Disposal container, Multi-Safe twin plus REC, 3,000 ml, grey, with plastic label, label/print: yellow/red, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 233 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 62 piece(s)/carton

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Absorbent liner	yes
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

Size

Max. volume	3000 ml
Diameter of base	112 mm
Height of product	233 mm

Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Disposal container, Multi-Safe twin plus REC,
3,000 ml

Packaging

Minimum order qty.	62
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	62
Piece(s) / pallet	744
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.241 kg
Weight of case	16.53 kg
EAN of inner box	4066916004002
EAN case	4066916004019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

