



Product Specification
Disposal container, Multi-Safe twin plus REC,
5,000 ml

Page 1

Product description

Order number 77.3884.250
Product description Disposal container, Multi-Safe twin plus REC, 5,000 ml, grey, with plastic label, label/print: yellow/red, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 379 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 34 piece(s)/carton

Product characteristics

Label/ Print with plastic label
Colour of print/label yellow/red
Example of use disposal of pointed and sharp objects
Absorbent liner yes
Shape of opening oval
Removal aid - turning aid Removal and rotation aid

Size

Max. volume 5000 ml
Diameter of base 112 mm
Height of product 379 mm

Material & colours

Product material recycled Polypropylene (recycled PP)
Colour of product grey
Closure material Polypropylene (PP)
Colour of cap red
Colour of lid red

Purity & certification

Satisfies the requirement TRBA 250, UN 3291, DIN EN ISO 23907
Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Disposal container, Multi-Safe twin plus REC,
5,000 ml

Packaging

Minimum order qty.	34
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	34
Piece(s) / pallet	408
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.342 kg
Weight of case	13.11 kg
EAN of inner box	4066916004125
EAN case	4066916004132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

