



Product description

| | |
|---------------------|---|
| Order number | 77.3884.350 |
| Product description | Disposal container, Multi-Safe twin plus REC, 5,000 ml, grey, with plastic label, label/print: yellow/red, biohazard labelling, material: recycled PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 379 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 34 piece(s)/carton |

Product characteristics

| | |
|---------------------------|---------------------------------------|
| Label/ Print | with plastic label |
| Colour of print/label | yellow/red |
| Example of use | disposal of pointed and sharp objects |
| Absorbent liner | yes |
| Shape of opening | oval |
| Removal aid - turning aid | Removal and rotation aid |

Size

| | |
|-------------------|---------|
| Max. volume | 5000 ml |
| Diameter of base | 112 mm |
| Height of product | 379 mm |

Material & colours

| | |
|-------------------|--------------------------------------|
| Product material | recycled Polypropylene (recycled PP) |
| Colour of product | grey |
| Closure material | Polypropylene (PP) |
| Colour of cap | red |
| Colour of lid | red |

Purity & certification

| | |
|---------------------------|-------------------------------------|
| Satisfies the requirement | TRBA 250, UN 3291, DIN EN ISO 23907 |
| Product category | no medical device IVD |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal container, Multi-Safe twin plus REC,
5,000 ml, biohazard labelling

Page 2

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 34 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 34 |
| Piece(s) / pallet | 408 |
| Depth of case | 800 mm |
| Width of case | 400 mm |
| Height of case | 380 mm |
| Case volume | 0.1216 cbm |
| Weight of product | 0.342 kg |
| Weight of case | 13.11 kg |
| EAN of inner box | 4066916004446 |
| EAN case | 4066916004453 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT