



## Product description

Order number	77.3890.202
Product description	Disposal container, Multi-Safe sani REC, 200 ml, grey, material: recycled PP, opening: oval, height: 131 mm, lid (LxW): 75 x 47 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 110 piece(s)/carton

## Product characteristics

Shape of opening	oval
Removal aid - turning aid	Removal aid

## Size

Max. volume	200 ml
Diameter of base	57 mm x 32 mm
Diameter of opening	33 mm x 23 mm
Height of product	131 mm

## Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Colour of cap	red
Colour of lid	red

## Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order qty.	110
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	110
Piece(s) / pallet	3300
Depth of case	397 mm
Width of case	397 mm
Height of case	325 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal container, Multi-Safe sani REC, 200 ml**

**Page 2**

Case volume	0.0512 cbm
Weight of product	0.0483 kg
Weight of case	6 kg
EAN of inner box	4066916004484
EAN case	4066916004491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**