



Product description

| | |
|---------------------|--|
| Order number | 77.3897.260 |
| Product description | Disposal container, Multi-Safe euroMatic® , 6,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labeling, material: PP, opening: oval, height: 330 mm, flat cap, lid (LxW): 164 x 164 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 40 piece(s)/carton |
| UN identification | yes |

Product characteristics

| | |
|---------------------------|--------------------|
| Label/ Print | with plastic label |
| Colour of print/label | yellow/red |
| Cap | flat cap |
| Shape of opening | oval |
| Removal aid - turning aid | Removal aid |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|---------------------|-----------------|
| Max. volume | 6000 ml |
| Diameter of base | 132 mm x 132 mm |
| Diameter of opening | 110 mm x 60 mm |
| Height of product | 330 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | yellow |
| Colour of cap | red |
| Colour of lid | red |

Purity & certification

| | |
|---------------------------|-------------------------------------|
| Satisfies the requirement | TRBA 250, UN 3291, DIN EN ISO 23907 |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 40 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 40 |
| Piece (s) / outer case | 40 |
| Piece(s) / pallet | 360 |
| Depth of case | 392 mm |
| Width of case | 514 mm |
| Height of case | 802 mm |
| Case volume | 0.1616 cbm |
| Weight of product | 0.4398 kg |
| Weight of case | 19 kg |
| EAN of inner box | 4038917962995 |
| EAN case | 4038917963008 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

