



### **Product description**

Order number	77.3897.260
Product description	Disposal container, Multi-Safe euroMatic® , 6,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labelling, material: PP, opening: oval, height: 330 mm, flat cap, lid (LxW): 164 x 164 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 40 piece(s)/carton
UN identification	yes

### **Product characteristics**

Label/ Print	with plastic label
Colour of print/label	yellow/red
Cap	flat cap
Shape of opening	oval
Removal aid - turning aid	Removal aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Size**

Max. volume	6000 ml
Diameter of base	132 mm x 132 mm
Diameter of opening	110 mm x 60 mm
Height of product	330 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	40
Type of smallest subpackaging	case
Piece(s) / inner box	40
Piece (s) / outer case	40
Piece(s) / pallet	360
Depth of case	392 mm
Width of case	514 mm
Height of case	802 mm
Case volume	0.1616 cbm
Weight of product	0.4398 kg
Weight of case	19 kg
EAN of inner box	4038917962995
EAN case	4038917963008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

