



### **Product description**

Order number	77.3898.106
Product description	Disposal container, Multi-Safe medi, 6 l, yellow, round, material: PP, with assembled bracket, opening: round, height: 195 mm, flat cap, Ø of cap: 230 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 480 piece(s)/carton
UN identification	yes

### **Product characteristics**

Cap	flat cap
Design	with assembled bracket
Shape of opening	round
Removal aid - turning aid	Removal and rotation aid

### **Size**

Max. volume	6 l
Diameter of base	190 mm
Diameter of lid	230 mm
Diameter of opening	88 mm
Height of product	195 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	480
Type of smallest subpackaging	pallet
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	480
Depth of case	792 mm
Width of case	1090 mm
Height of case	1192 mm
Case volume	1.029 cbm
Weight of product	0.3126 kg
Weight of case	160 kg
EAN of inner box	4038917543422
EAN case	4038917543439

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

