



### **Product description**

Order number	51.1534.100
Product description	Salivette® VD, with cotton swab, cap: light blue, with paper label, label/print: light blue/white, for Saliva examination for virus diagnostics, 100 piece(s)/bag

### **Product characteristics**

Label/ Print	with paper label
Colour of print/label	light blue/white
Closure type	push cap
Design	with cotton swab
Application	Saliva examination for virus diagnostics
Recommended centrifugation	500 x g
Centrifugation time	1 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	16.8 mm
Length including cap	97 mm
Length excluding cap	84 mm

### **Material & colours**

Product material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	light blue

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	28000
Depth of box	226 mm
Width of box	101 mm
Height of box	119 mm
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0069 kg
Weight of case	4.09 kg
EAN of inner box	4038917581882
EAN case	4038917581875

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

