



### **Product description**

Order number	95.1123
Product description	Cold transport container, english

### **Product characteristics**

Label/ Print	with print
Colour of print/label	red
Application	Cold shipping
Absorbent liner	no

### **Size**

Width of product	71 mm
Height of product	170 mm
Length of product	140 mm

### **Material & colours**

Product material	Acrylnitrile butadiene styrene (ABS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	6
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	6
Piece(s) / pallet	648
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.3448 kg
Weight of case	2.3 kg
EAN of inner box	4038917085403
EAN case	4038917004084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

