



Product description

Order number 95.995
Product description Cold transport container, german

Product characteristics

Label/ Print with print
Colour of print/label red
Application Cold shipping
Absorbent liner no

Size

Width of product 71 mm
Height of product 170 mm
Length of product 140 mm

Material & colours

Product material Acrylnitrile butadiene styrene (ABS)
Colour of product white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	6
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	6
Piece(s) / pallet	648
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.3448 kg
Weight of case	2.3 kg
EAN of inner box	4038917084901
EAN case	4038917003667

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

