



## Product Specification

Transport packaging, LabMailer Small, with absorbent liner, length: 82 mm, Ø opening: 78 mm, cap assembled

Page 1



### Product description

Order number	78.9925.418
Product description	Transport packaging, LabMailer Small, with absorbent liner, (LxW): 82 x 82 mm, Ø opening: 78 mm, cap assembled, cap: red

### Product characteristics

Application	Transportation of diagnostic specimens
Absorbent liner	yes

### Size

Diameter	82 mm
Diameter of opening	78 mm
Width of product	82 mm
Length of product	82 mm
Length including cap	179 mm
Usable length	153 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Transport packaging, LabMailer Small, with  
absorbent liner, length: 82 mm, Ø opening: 78  
mm, cap assembled

Page 2

### Material & colours

Product material	cardboard
Colour of cap	red

### Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	50
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	1600
Case volume	0.0692 cbm
Weight of case	11 kg
EAN of inner box	4038917155526
EAN case	4038917039437

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**