



Product description

Order number	86.974
Product description	Mailing envelope, material: PE, transparent, (LxW): 440 x 320 mm, 300 piece(s)/case

Product characteristics

Absorbent liner	no
-----------------	----

Size

Width of product	320 mm
Length of product	440 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 300
Type of smallest case
subpackaging
Piece(s) / inner box 300
Piece (s) / outer case 300
Piece(s) / pallet 7200
Depth of case 295 mm
Width of case 320 mm
Height of case 395 mm
Case volume 0.0373 cbm
Weight of product 0.045 kg
Weight of case 14.2 kg
EAN of inner box 4038917084758
EAN case 4038917003537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

