



Product description

Order number	95.978
Product description	Mailing envelope, material: PE, (LxW): 230 x 300 mm, with 50 mm adhesive tape, film thickness: 0.045 mm, 100 piece(s)/case

Product characteristics

Design	with 50 mm adhesive tape
Absorbent liner	no

Size

Width of product	300 mm
Length of product	230 mm

Material & colours

Product material	Polyethylene (PE)
------------------	-------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging case
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 56000
Depth of case 290 mm
Width of case 145 mm
Height of case 400 mm
Case volume 0.0168 cbm
Weight of product 0.0057 kg
Weight of case 6.12 kg
EAN of inner box 4038917224178
EAN case 4038917224161

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

