



### **Product description**

Order number	95.1712
Product description	Superabsorber pack, suitable for transport case T 15 and B 17, (LxW): 400 x 300 mm, fluid intake: 1.2 l, 100 piece(s)/case

### **Product characteristics**

Application	Mailing of samples, Multi-Safe disposal boxes
Suitable for	transport case T 15 and B 17

### **Size**

Width of product	300 mm
Length of product	400 mm

### **Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Superabsorber pack, suitable for transport case**  
**T 15 and B 17, (LxW): 400 x 300 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      100  
Type of smallest  
subpackaging                              case  
Piece(s) / inner box                      100  
Piece (s) / outer case                      100  
Piece(s) / pallet                            72000  
Depth of case                              420 mm  
Width of case                                320 mm  
Height of case                              140 mm  
Case volume                                0.0188 cbm  
Weight of product                        0.0283 kg  
Weight of case                              3.5 kg  
EAN of inner box                          4038917415774  
EAN case                                      4038917065894

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

