



### **Product description**

Order number 95.1717  
Product description Transport case B 17, with lateral rim for block rack,  
(LxWxH): 395 x 395 x 435 mm

### **Product characteristics**

Design with lateral rim for block rack

### **Size**

Width of product 395 mm  
Height of product 435 mm  
Length of product 395 mm

### **Material & colours**

Product material High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	18
Depth of case	403 mm
Width of case	403 mm
Height of case	460 mm
Case volume	0.0747 cbm
Weight of product	3.561 kg
Weight of case	4.5 kg
EAN of inner box	4038917087544
EAN case	4038917065610

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

