



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 95.1717.011   |
| Product description | Transport case B 13, flat design for block rack, (LxWxH):<br>395 x 395 x 300 mm |

### **Product characteristics**

|        |                            |
|--------|----------------------------|
| Design | flat design for block rack |
|--------|----------------------------|

### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 395 mm |
| Height of product | 300 mm |
| Length of product | 395 mm |

### **Material & colours**

|                  |                                   |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
|------------------|-----------------------------------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR              |
| Product category          | no medical device   IVD |
| Batched                   | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 1             |
| Piece(s) / pallet             | 30            |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 340 mm        |
| Case volume                   | 0.0539 cbm    |
| Weight of product             | 2.763 kg      |
| Weight of case                | 3.3 kg        |
| EAN of inner box              | 4038917133197 |
| EAN case                      | 4038917069250 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

