



Product description

| | |
|---------------------|---|
| Order number | 95.1720 |
| Product description | Transport case I 100, internal case with attachment, (LxWxH): 360 x 310 x 275 mm |

Product characteristics

| | |
|---------------------|-------------------------------|
| Design | internal case with attachment |
| Shape of inner tube | square |

Size

| | |
|-------------------|--------|
| Width of product | 310 mm |
| Height of product | 275 mm |
| Length of product | 360 mm |

Material & colours

| | |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
|------------------|-----------------------------------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | TRBA 100 |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 30 |
| Depth of case | 398 mm |
| Width of case | 358 mm |
| Height of case | 325 mm |
| Case volume | 0.0463 cbm |
| Weight of product | 1.613 kg |
| Weight of case | 2.05 kg |
| EAN of inner box | 4038917087551 |
| EAN case | 4038917076210 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

