




Product Specification
Screw cap, yellow, suitable for Mailing container
126 x 30 mm

Page 1



Product description

Order number 65.679.004
Product description Screw cap, yellow, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag

Product characteristics

Application Mailing of samples
Suitable for Mailing container 126 x 30 mm

Size

Length of product 20 mm

Material & colours

Product material Polypropylene (PP)
Colour of product yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap, yellow, suitable for Mailing container
126 x 30 mm

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 250
Type of smallest
subpackaging bag
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 41250
Depth of case 338 mm
Width of case 228 mm
Height of case 120 mm
Case volume 0.0092 cbm
Weight of product 0.0035 kg
Weight of case 1.1 kg
EAN of inner box 4038917115834
EAN case 4038917014939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT