



### **Product description**

Order number	78.670
Product description	Absorbent liner, suitable for Mailing container 126 x 30 mm, (LxW): 75 x 90 mm, 1,000 piece(s)/case

### **Product characteristics**

Application	Mailing of samples, Multi-Safe disposal boxes
Suitable for	Mailing container 126 x 30 mm

### **Size**

Width of product	90 mm
Length of product	75 mm

### **Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for Mailing container 126**  
**x 30 mm, (LxW): 75 x 90 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      4000  
Type of smallest  
subpackaging                              case  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      4000  
Piece(s) / pallet                            72000  
Depth of case                              515 mm  
Width of case                                360 mm  
Height of case                              360 mm  
Case volume                                0.0667 cbm  
Weight of product                        0.001 kg  
Weight of case                              4.36 kg  
EAN of inner box                         4038917377492  
EAN case                                      4038917159937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

