



Product description

Order number	78.671
Product description	Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm, 100 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposal boxes
Suitable for	cold transport container

Size

Width of product	55 mm
Length of product	80 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2500
Type of smallest subpackaging case
Piece(s) / inner box 100
Piece (s) / outer case 2500
Piece(s) / pallet 270000
Depth of case 383 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.014 cbm
Weight of product 0.0009 kg
Weight of case 2.57 kg
EAN of inner box 4038917083263
EAN case 4038917079860

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

