



Product Specification

Mailing container, transparent, construction: round, with absorbent liner, length: 126 mm, Ø opening: 30 mm, without cap

Page 1



Product description

Order number	78.898
Product description	Mailing container, round, transparent, construction: round, with absorbent liner, length: 126 mm, Ø opening: 30 mm, without cap, 50 piece(s)/bag

Product characteristics

Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	yes
Shape of opening	round
Liquid-tight	yes

Size

Diameter of opening	30 mm
Length of product	126 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Mailing container, transparent, construction:
round, with absorbent liner, length: 126 mm, Ø
opening: 30 mm, without cap**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement P 650, ADR
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 250
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 7500
Depth of case 398 mm
Width of case 358 mm
Height of case 325 mm
Case volume 0.0463 cbm
Weight of product 0.0145 kg
Weight of case 4.2 kg
EAN of inner box 4038917084543
EAN case 4038917003353

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT