



Product Specification
UriSet 24, with stabiliser, 3 l, with
Urine-Monovette® 10 ml, brown, with light
protection, graduated



Product description

| | |
|---------------------|---|
| Order number | 77.578.252 |
| Product description | UriSet 24, Complete set for optimum 24-hour urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, with stabiliser, brown, with light protection, graduated, material: PE, cap: green, with Urine-Monovette® 10 ml, 30 piece(s)/case |

Product characteristics

| | |
|---------------------|---|
| Type of preparation | with stabiliser |
| Graduation | yes |
| Closure type | screw cap |
| Light protection | yes |
| Design | with Urine-Monovette® 10 ml |
| Application | Complete set for optimum 24-hour urine collection |
| Viewing window | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
UriSet 24, with stabiliser, 3 l, with
Urine-Monovette® 10 ml, brown, with light
protection, graduated

Page 2

Size

| | |
|----------------------|-------------|
| Sample volume | 3.000.00 ml |
| Max. volume | 3 l |
| Diameter of opening | 83,5 mm |
| Width of product | 130 mm |
| Length of product | 125 mm |
| Length including cap | 246 mm |
| Length excluding cap | 243 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polyethylene (PE) |
| Colour of product | brown |
| Closure material | Polypropylene (PP) |
| Colour of cap | green |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 30 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 30 |
| Piece (s) / outer case | 30 |
| Piece(s) / pallet | 270 |
| Depth of case | 643 mm |
| Width of case | 403 mm |
| Height of case | 515 mm |
| Case volume | 0.1335 cbm |
| Weight of product | 0.2296 kg |
| Weight of case | 7.69 kg |
| EAN of inner box | 4038917113878 |
| EAN case | 4038917072762 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT