



Product description

Order number	75.1337.505
Product description	Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, Ø opening: 37 mm, green, graduated, material: PP, 75 piece(s)/bag

Product characteristics

Graduation	yes
Volume of work	30 ml

Size

Sample volume	30 ml
Max. volume	30 ml
Diameter	37 mm
Diameter of opening	37 mm
Height of product	40 mm
Length of product	40 mm
Length excluding cap	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, graduated, PP, green

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product green

Purity & certification

Product category Medical device
CE certified CE supplier
Batched yes

Packaging

Minimum order qty. 3750
Type of smallest subpackaging bag
Piece(s) / inner box 3750
Piece (s) / outer case 3750
Piece(s) / pallet 180000
Depth of case 190 mm
Width of case 390 mm
Height of case 390 mm
Case volume 0.0289 cbm
Weight of product 0.0015 kg
Weight of case 6 kg
EAN of inner box 4038917127752
EAN case 4042836002069

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT