



Product description

| | |
|---------------------|--|
| Order number | 75.1355 |
| Product description | Urine container, collection and storage of urine, max. working volume: 125 ml, (ØxH): 66 x 67 mm, Ø opening: 66 mm, highly transparent, graduated, material: PP, 50 piece(s)/bag |

Product characteristics

| | |
|----------------|---------------------------------|
| Graduation | yes |
| Application | collection and storage of urine |
| Volume of work | 125 ml |

Size

| | |
|----------------------|--------|
| Sample volume | 125 ml |
| Max. volume | 125 ml |
| Diameter | 66 mm |
| Diameter of opening | 66 mm |
| Height of product | 67 mm |
| Length excluding cap | 67 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 125 ml, (ØxH): 66 x 67 mm, PP,
highly transparent**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product highly transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE supplier
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 1000
Piece(s) / pallet 24000
Depth of case 390 mm
Width of case 265 mm
Height of case 590 mm
Case volume 0.061 cbm
Weight of product 0.0057 kg
Weight of case 6.2 kg
EAN of inner box 4038917409568
EAN case 4042836010088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT