



### **Product description**

Order number	75.560.901
Product description	Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS

### **Product characteristics**

Graduation	no
Volume of work	148 ml

### **Size**

Sample volume	148 ml
Diameter	73 mm
Height of product	67 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                    yes

**Packaging**

Minimum order qty.                    2500  
Piece(s) / inner box                    100  
Piece (s) / outer case                   2500  
Piece(s) / pallet                        90000  
Depth of case                            457 mm  
Width of case                            368 mm  
Height of case                           356 mm  
Case volume                              0.0599 cbm  
Weight of product                       0.0028 kg  
Weight of case                           8.16 kg  
EAN of inner box                        4038917318198  
EAN case                                  4038917916547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

