



Product Specification

Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label

Page 1



Product description

Order number	75.9922.721
Product description	Multi-purpose container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	70 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label

Page 2

Size

Sample volume	70 ml
Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm
Length of product	55 mm
Length excluding cap	55 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	blister
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	3200
Depth of case	548 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0834 cbm
Weight of product	0.0112 kg
Weight of case	3.5 kg
EAN of inner box	4038917152525
EAN case	4038917347181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT