



Product description

Order number	76.1340.500
Product description	Lid, for 30-ml medicine cup, material: PE, white, 750 piece(s)/bag

Product characteristics

Example of use	Oral administration of medications in liquid or tablet form
Suitable for	30-ml medicine cup

Material & colours

Product material	Polyethylene (PE)
Colour of product	white
Colour of cap	white

Purity & certification

Product category	Medical device
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	180000
Depth of case	190 mm
Width of case	390 mm
Height of case	390 mm
Case volume	0.0289 cbm
Weight of product	0.0007 kg
Weight of case	2.97 kg
EAN of inner box	4038917127790
EAN case	4042836007149

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

