



Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile



Product description

Order number	51.9923.820
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, skirted conical base, transparent, cap assembled, with paper label, preparation: Boric acid, non-sterile, 500 piece(s)/case

Product characteristics

Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile

Page 2

Size

Sample volume	25 ml
Diameter	25 mm
Length of product	90 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	red

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

Packaging

Minimum order qty.	10000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	15000
Case volume	0.059 cbm
Weight of product	0.009 kg
Weight of case	5.18 kg
EAN of inner box	4038917154550
EAN case	4038917030090

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT