



### **Product description**

Order number 60.9921.830  
Product description Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PS, skirted conical base, transparent, screw cap, yellow, cap assembled, 100 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type screw cap  
Base shape skirted conical base  
Volume of work 10 ml  
Centrifugation max (RCF) 4000 x g

### **Size**

Diameter 16 mm  
Length of product 97 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### **Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0079 kg
Weight of case	8.7 kg
EAN of inner box	4038917151139
EAN case	4038917027717

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

