



Product description

Order number	10.250.001
Product description	Urine Monovette® , nominal volume: 3.2 ml, without preparation, Urine tube, (LxØ): 75 x 13 mm, transparent, fill mark, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	3,2 ml
Recommended centrifugation	400 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 3.2 ml, cap yellow, (LxØ): 75 x 13 mm, 64 piece(s)/bag

Page 2

Centrifugation time 5 Min
Connection Luer (m)

Size

Sample volume 3.2 ml
Diameter 13 mm
Length including cap 90 mm
Length excluding cap 75 mm
Length including cap and plunger 109 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 512
Type of smallest subpackaging bag
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 36864
Depth of case 388 mm
Width of case 358 mm
Height of case 150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Urine Monovette® , 3.2 ml, cap yellow, (LxØ): 75 x 13 mm, 64 piece(s)/bag

Page 3

Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.23 kg
EAN of inner box	4038917279499
EAN case	4038917279482

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT