



Product description

| | |
|---------------------|---|
| Order number | 10.252.030 |
| Product description | Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, white, with light protection, graduated, material: PP, Luer screw cap, cap yellow, with plastic label, 64 piece(s)/bag |

Product characteristics

| | |
|-----------------------|--|
| Type of collection | Aspiration |
| Type of preparation | without preparation |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/yellow |
| Graduation | yes |
| Cap | Luer screw cap |
| Closure type | screw cap |
| Light protection | yes |
| Application | Urine tube |
| Example of use | Urine collection, transport, and storage |
| Base shape | round base |
| Volume of work | 10 ml |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 2

Recommended centrifugation 400 x g
Centrifugation time 5 Min
Connection Luer (m)

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 117 mm
Length excluding cap 102 mm
Length including cap and plunger 127 mm

Material & colours

Product material Polypropylene (PP)
Colour of product white
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 512
Type of smallest subpackaging bag
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 24576
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 3

| | |
|-------------------|---------------|
| Case volume | 0.0315 cbm |
| Weight of product | 0.0085 kg |
| Weight of case | 4.85 kg |
| EAN of inner box | 4038917279987 |
| EAN case | 4038917279840 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT