



Product description

Order number	10.252.388
Product description	Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix, sterile, 1 piece(s)/blister

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 2

Volume of work	10 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

Size

Sample volume	10 ml
Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 3

Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0084 kg
Weight of case	5.82 kg
EAN of inner box	4038917949828
EAN case	4038917949835

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT