



Product Specification

Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 1



Product description

Order number	10.253.001
Product description	Urine Monovette® , Boric acid, nominal volume: 10 ml, Boric acid, Urine tube, (LxØ): 102 x 15 mm, transparent, material: PP, Luer screw cap, cap green, incl. collection tip, with paper label, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	no
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag

Connection Luer (m)

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 117 mm
Length excluding cap 102 mm
Length including cap and plunger 127 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Safety data sheet yes
Batched yes

Packaging

Minimum order qty. 512
Type of smallest subpackaging bag
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 24576
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 3

Weight of product	0.0086 kg
Weight of case	4.9 kg
EAN of inner box	4038917280013
EAN case	4038917279871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT