



## Product Specification

Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 1



### Product description

Order number	10.253.021
Product description	Urine Monovette® , Boric acid, nominal volume: 10 ml, Boric acid, Urine tube, (LxØ): 102 x 15 mm, transparent, material: PP, Luer screw cap, cap green, with paper label, sterile, 1 piece(s)/blister

### Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	no
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 2

Connection Luer (m)

### Size

Sample volume 10 ml  
Diameter 15 mm  
Length including cap 117 mm  
Length excluding cap 102 mm  
Length including cap and plunger 127 mm

### Material & colours

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap green  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Safety data sheet yes  
Batched yes

### Packaging

Minimum order qty. 20000  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 500  
Piece(s) / pallet 13500  
Depth of box 351 mm  
Width of box 151 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

**Page 3**

Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0085 kg
Weight of case	5.87 kg
EAN of inner box	4038917323864
EAN case	4038917323857

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**