



## Product Specification

**Urine Monovette® , Boric acid, 3.2 ml, cap green,  
(LxØ): 75 x 13 mm, 1 piece(s)/blister**

**Page 1**

### Product description

Order number	10.256.021
Product description	Urine Monovette® , Boric acid, nominal volume: 3.2 ml, Boric acid, Urine tube, (LxØ): 75 x 13 mm, transparent, fill mark, material: PP, Luer screw cap, cap green, incl. collection tip, with paper label, suitable for Tempus600® sample transport system, sterile, 1 piece(s)/blister

### Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus600® sample transport system
Base shape	round base
Volume of work	3,2 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

### Size

Sample volume	3.2 ml
Diameter	13 mm
Length including cap	90 mm
Length excluding cap	75 mm
Length including cap and plunger	109 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine Monovette® , Boric acid, 3.2 ml, cap green,  
(LxØ): 75 x 13 mm, 1 piece(s)/blister**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Safety data sheet	yes
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0046 kg
Weight of case	4.02 kg
EAN of inner box	4038917983143
EAN case	4038917983150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**