



### Product description

Order number	11.2252.001
Product description	V-Monovette® Urine, nominal volume: 4 ml, without preparation, vacuum system for closed urine transfer, (LxØ): 75 x 13 mm, transparent, fill mark, material: PET, membrane cap, cap yellow, with paper label, sterile, 50 piece(s)/bag

### Product characteristics

Type of collection	Vacuum
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	round base
Volume of work	4 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**V-Monovette® Urine, 4 ml, cap yellow, (LxØ): 75 x 13 mm, 50 piece(s)/bag**

**Page 2**

Recommended centrifugation 400 x g  
Centrifugation time 5 Min  
Connection V-Monovette®

**Size**

Sample volume 4 ml  
Diameter 13 mm  
Length including cap 81 mm  
Length excluding cap 75 mm

**Material & colours**

Product material Polyethylene therephthalate (PET)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Membrane material rubber  
Tube material Polyethylene therephthalate (PET)

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest subpackaging bag  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of case 398 mm  
Width of case 368 mm  
Height of case 130 mm  
Case volume 0.019 cbm  
Weight of product 0.006 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**V-Monovette® Urine, 4 ml, cap yellow, (LxØ): 75 x 13 mm, 50 piece(s)/bag**

Weight of case	3.46 kg
EAN of inner box	4038917244824
EAN case	4038917431842

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

