



Product description

Order number	11.2253.001
Product description	V-Monovette® Urine, Boric acid, nominal volume: 4 ml, Boric acid, vacuum system for closed urine transfer, (LxØ): 75 x 13 mm, transparent, fill mark, material: PET, membrane cap, cap green, with paper label, sterile, 50 piece(s)/bag

Product characteristics

Type of collection	Vacuum
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**V-Monovette® Urine, Boric acid, 4 ml, cap green,
(LxØ): 75 x 13 mm, 50 piece(s)/bag**

Page 2

Volume of work 4 ml
Connection V-Monovette®

Size

Sample volume 4 ml
Diameter 13 mm
Length including cap 81 mm
Length excluding cap 75 mm

Material & colours

Product material Polyethylene terephthalate (PET)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Membrane material rubber
Tube material Polyethylene terephthalate (PET)

Purity & certification

Product category In vitro diagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of case 398 mm
Width of case 368 mm
Height of case 130 mm
Case volume 0.019 cbm
Weight of product 0.006 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
V-Monovette® Urine, Boric acid, 4 ml, cap green,
(LxØ): 75 x 13 mm, 50 piece(s)/bag

Page 3

Weight of case	3.47 kg
EAN of inner box	4038917245128
EAN case	4038917245111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT