



Product description

Order number	14.1513
Product description	Colour-coded cap, for S-Monovette [®] , pink, 100 piece(s)/bag

Product characteristics

Suitable for	S-Monovette [®]
Connection	S-Monovette [®]

Size

Length of product	13 mm
-------------------	-------

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	pink

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 324000
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.0003 kg
Weight of case 0.47 kg
EAN of inner box 4038917086790
EAN case 4038917073943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

