



Product Specification
Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation



Product description

| | |
|---------------------|--|
| Order number | 53.419 |
| Product description | Seraplas® valve filter, V13, blue, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case |

Product characteristics

| | |
|--------------|--|
| Application | venous blood collection |
| Suitable for | separation of serum/plasma from the blood cells after centrifugation |

Size

| | |
|-------------------|-------|
| Length of product | 92 mm |
|-------------------|-------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation

Material & colours

Product material Polyethylene (PE)
Colour of product blue

Purity & certification

Product category no medical device | IVD
Safety data sheet no
Batched no

Packaging

Minimum order qty. 1000
Type of smallest subpackaging case
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 75000
Depth of box 119 mm
Width of box 99 mm
Height of box 98 mm
Depth of case 513 mm
Width of case 293 mm
Height of case 140 mm
Case volume 0.021 cbm
Weight of product 0.0016 kg
Weight of case 2.15 kg
EAN of inner box 4038917081481
EAN case 4038917000574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

