



Product Specification
Seraplas® valve filter, green, for separation of serum/plasma from the blood cells after centrifugation

Page 1



Product description

Order number	53.420
Product description	Seraplas® valve filter, V16, green, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

Product characteristics

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

Size

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Seraplas® valve filter, green, for separation of serum/plasma from the blood cells after centrifugation

Page 2

Material & colours

Product material Polyethylene (PE)
Colour of product green

Purity & certification

Product category no medical device | IVD
Safety data sheet no
Batched no

Packaging

Minimum order qty. 1000
Type of smallest subpackaging case
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 54000
Depth of box 156 mm
Width of box 99 mm
Height of box 130 mm
Depth of case 513 mm
Width of case 328 mm
Height of case 150 mm
Case volume 0.0252 cbm
Weight of product 0.0023 kg
Weight of case 2.98 kg
EAN of inner box 4038917081498
EAN case 4038917000581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT