



### Product description

Order number	55.1570.002
Product description	Adapter tube, (LxØ): 54 x 11 mm, material: PP, brown, 250 piece(s)/bag

### Product characteristics

Light protection	yes
Suitable for	S-Monovette® with Ø 8 mm
Base shape	flat base

### Size

Diameter	11 mm
Length of product	54 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      1000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      250  
Piece (s) / outer case                      1000  
Piece(s) / pallet                            78000  
Depth of case                              428 mm  
Width of case                                323 mm  
Height of case                               130 mm  
Case volume                                0.018 cbm  
Weight of product                         0.0012 kg  
Weight of case                              1.54 kg  
EAN of inner box                          4038917131018  
EAN case                                    4038917023573

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

