



## Product Specification

Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
cap assembled, 100 piece(s)/bag

Page 1



### Product description

Order number	62.617
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Graduation	yes
Base shape	conical false bottom, rounded tube bottom
Volume of work	3,5 ml

### Size

Diameter	13 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
cap assembled, 100 piece(s)/bag**

**Page 2**

### Material & colours

Product material Polypropylene (PP)  
Colour of product transparent

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### Packaging

Minimum order qty. 1000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 36000  
Depth of case 458 mm  
Width of case 408 mm  
Height of case 205 mm  
Case volume 0.0383 cbm  
Weight of product 0.0036 kg  
Weight of case 4.34 kg  
EAN of inner box 4038917082921  
EAN case 4038917079150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**